

Installation Guide for ByoDyn version 4.8

Adrián López García de Lomana, Alex Gómez-Garrido, Miguel Hernández, Pau Rué
Queralt and Jordi Villà-Freixa

— December 17, 2008 —

This document specifies the requirements and different steps to install
ByoDyn on a Linux box (Fedora Core 2, 4 and 6).

Computational Biochemistry and Biophysics Laboratory
Research Unit on Biomedical Informatics
Universitat Pompeu Fabra - Institut Municipal d'Investigació Mèdica
c/ Dr. Aiguader 88, 08003, Barcelona, Spain
<http://cbbl.imim.es>

Contents

1	ByoDyn Requirements	3
1.1	Mandatory	3
1.2	Optional	4
2	Installation of ByoDyn	5

1 ByoDyn Requirements

1.1 Mandatory

- **Python:** version 2.5 or newer is required <http://www.python.org/download>.
- **SciPy:** version 0.6.0 or newer available at <http://www.scipy.org/>. To install SciPy you also need **numpy** available at <http://numeric.scipy.org/>. Please install **numpy** version 1.1.0 or newer. For an smooth use of SciPy, make sure you have built BLAS and LAPACK libraries as explained at http://www.scipy.org/Installing_SciPy/BuildingGeneral.
- **libSBML:** version 3.2.0 or newer available at <http://www.sbml.org/libsbml.html>. Make sure you use at least the Python binder adding the flag `--with-python` while configuring.
- **PORT Library:** from Netlib. Please follow the required specified steps for the correct functioning of the local search routines of ByoDyn:

- Use the following command to download the code:

```
rsync -avz netlib.org::netlib/port .
```

- Modify the Makefile:

- * If your Fortran compiler is not `f77` add the following line at the top the file:

```
FC=gfortran
```

`gfortran` is the default Fortran compiler in the latest Fedora distribution. If you have another compiler, change `gfortran` by the one you have.

- * Change the line

```
LIB=portP
```

by

```
LIB=port
```

we are simply changing the name of the library.

- * Remove the line

```
n5err.o\
```

the file `n5err.f` was not found in the server and is not required in this case

- * Change the line

```
update lib$(LIB).a $?
```

by

```
ar rcs lib$(LIB).a $?
```

- * Change the line

```
FFLAGS=-O
```

by

`FFLAGS=-O -fPIC`

if your architecture is 64 bits.

- Execute `make`.

Once you compile the source code, the file called `libport.a` has to be copied to the `byodyn/lib/local_search/PORT` directory.

- **Gnuplot**: Gnuplot 3.7 or newer needed to create the output graphs. Available at <http://www.gnuplot.i>
- **matplotlib**: Version 0.91.1 or newer is required. You can download it from [http://matplotlib.sourcef](http://matplotlib.sourceforge). It is required for the plots output by ByoDyn.

1.2 Optional

In addition to this, other pieces of software add extra functionality to the ByoDyn platform:

- **Doxygen**: required to build the API of ByoDyn in HTML and PDF formats. It has been tested under Doxygen version 1.5.2.
- **Octave**: is an optional tool to substitute SciPy functions for the integration of systems of ordinary differential equations (ODEs). On the contrary Octave is required for the solution of differential-algebraic equations (DAEs). Moreover Octave is compulsory required for clustering. Octave version 3.0.1 or newer is required. Octave is available at <http://www.octave.org/>.
- **OpenModelica**: version 1.4.3 or newer is required for the resolution of ODEs, DAEs and DAEs with events. Information for the installation can be found at <http://www.ida.liu.se/~pelab/mod>. Some environmental variables are required:
 - `OPENMODELICAHOME`: pointing to the `openmodelica` directory.
 - `PATH`: `omc` and `./` should be on the `PATH` for a correct execution of the software.
- **XPP/XPPAUT**: version 5.91 or newer is required for solving SBML models with delays. Information for the installation can be found at <http://www.math.pitt.edu/~bard/xpp/xpp.html>.
- **Open MPI**: the command `mpirun` is required to run ByoDyn in parallel. Version 1.2.3 or higher should be available.
- **ScientificPython**: the module `Scientific.MPI` is required for launching parallel calculations. Version 2.6 or newer is required. Please check <http://dirac.cnrs-orleans.fr/plone/software/> for further information. Be sure you also create and call the executable `mpipython`. To build it please follow the instructions described in the `README.MPI` file of Scientific. If you will run ByoDyn on a single processor machine, the installation of this software is not required.

2 Installation of ByoDyn

Once all those software pieces have been installed into the machine and are accessible by ByoDyn you should set up a few things before you run it.

- First, add the ByoDyn executable to your path. For bash this would be done by:
 - edit your `$HOME/.bashrc`.
 - define the ByoDyn path: `export BYODYN_PATH=$HOME/where_ByoDyn_is`
For example if you have ByoDyn at `$HOME/a_directory/another_directory` you should add the line
`export BYODYN_PATH=$HOME/a_directory/another_directory/byodyn`
at the `.bashrc` file.
 - Add ByoDyn executable file to the path: `export PATH="$BYODYN_PATH/bin:$PATH"`
- Make sure that the file `libport.a` is at the correct directory as specified at the Section [1.1](#).
- Finally, execute the command `byodyn`. The following lines should prompt:
`ByoDyn version 4.6 is running ...`
`Run byodyn -h for help.`
`... exiting from ByoDyn.`

At the case of a newly installed ByoDyn, we suggest to run the ByoDyn tests in order check the correct functioning of the program. For that, type `byodyn --testing` on your terminal.